

# STAR SERIES



Fire Resistant Polymer Body | ClearVision Technology



## Star 105+105

Intensity	1,80,000 + 1,80,000 Lux
Shadowless Feature	Yes
Numbers of LED's	105 + 105
Intensity Control	10~100 Digital Capacitive Touch (Fully Remote Controlled)
Fitting Options Available	Ceiling 180 Degree, Ceiling 360 Degree
Spot Diameter	100-150 mm
Depth of Light Penetration	150 mm
Colour Temperature	3800 – 4800 K
Colour Rendering Index	95 Ra
Focusing	Adjustable
LED Average life	>50,000 Hrs.
Diameter of Light	700 mm + 700 mm
Power Supply	220V/50Hz AC
Power Consumption	100 + 100 Watt

## Features

Shadowless Light, Memory Function, High quality sterilizable handle, With protective fire resistant Polycarbonate Body, fast disinfection process, protects UV and IR rays and improve light penetration. Aerodynamic design for ventilation of heat. Body is made of high quality fire resistant polymer



Memory Function



Endo Mode



Battery Backup: Optional

## Star 105

Intensity	1,80,000 Lux
Shadowless Feature	Yes
Numbers of LED's	105
Intensity Control	10~100 Digital Capacitive Touch (Fully Remote Controlled)
Fitting Options Available	Ceiling 180 Degree, Ceiling 360 Degree, Mobile Stand, Wall Mounted
Spot Diameter	100-150 mm
Depth of Light Penetration	150 mm
Colour Temperature	3800 – 4800 K
Colour Rendering Index	95 Ra
Focusing	Adjustable
LED Average life	>50,000 Hrs.
Diameter of Light	700 mm
Power Supply	220V/50Hz AC
Power Consumption	100 Watt

## Features

Shadowless Light, Memory Function, High quality sterilizable handle, With protective fire resistant Polycarbonate Body, fast disinfection process, protects UV and IR rays and improve light penetration. Aerodynamic design for ventilation of heat. Body is made of high quality fire resistant polymer



Memory Function



Endo Mode



Battery Backup: Optional

## Star 84+84

Intensity	1,40,000 + 1,40,000 Lux
Shadowless Feature	Yes
Numbers of LED's	84 + 84
Intensity Control	10~100 Digital Capacitive Touch (Fully Remote Controlled)
Fitting Options Available	Ceiling 180 Degree, Ceiling 360 Degree
Spot Diameter	100-150 mm
Depth of Light Penetration	150 mm
Colour Temperature	3800 – 4800 K
Colour Rendering Index	95 Ra
Focusing	Adjustable
LED Average life	>50,000 Hrs.
Diameter of Light	700 mm + 700 mm
Power Supply	220V/50Hz AC
Power Consumption	100 + 100 Watt

## Features

Shadowless Light, Memory Function, High quality sterilizable handle, With protective fire resistant Polycarbonate Body, fast disinfection process, protects UV and IR rays and improve light penetration. Aerodynamic design for ventilation of heat. Body is made of high quality fire resistant polymer



Memory Function



Endo Mode



Battery Backup: Optional

## Star 84

Intensity	1,40,000 Lux
Shadowless Feature	Yes
Numbers of LED's	84
Intensity Control	10~100 Digital Capacitive Touch (Fully Remote Controlled)
Fitting Options Available	Ceiling 180 Degree, Ceiling 360 Degree, Mobile Stand, Wall Mounted
Spot Diameter	100-150 mm
Depth of Light Penetration	150 mm
Colour Temperature	3800 – 4800 K
Colour Rendering Index	95 Ra
Focusing	Adjustable
LED Average life	>50,000 Hrs.
Diameter of Light	700 mm
Power Supply	220V/50Hz AC
Power Consumption	100 Watt

## Features

Shadowless Light, Memory Function, High quality sterilizable handle, With protective fire resistant Polycarbonate Body, fast disinfection process, protects UV and IR rays and improve light penetration. Aerodynamic design for ventilation of heat. Body is made of high quality fire resistant polymer



Memory Function



Endo Mode



Battery Backup: Optional

## Star 63 + 63

Intensity	1,20,000 + 1,20,000 Lux
Shadowless Feature	Yes
Numbers of LEDs	63 + 63 (Warm + Cool White)
Intensity Control	10~100 Digital Capacitive Touch (Fully Remote Controlled)
Fitting Options Available	Ceiling Mounted - 340 Degrees/ 360 Degrees
Spot Diameter	100-150 mm
Depth of Light Penetration	150 mm
Colour Temperature	3800 – 4800K
Colour Rendering Index	93 Ra
Focusing	Adjustable
LED Average life	>50,000 Hrs.
Diameter of Light	650 mm + 650 mm
Power Supply	220V/50Hz AC
Power Consumption	70 + 70 Watt
Optional Special Features	Mobile App Connectivity, External Wall Panel, Battery Backup

## Features

Shadowless Light, Memory Function, High quality sterilizable handle, With protective fire resistant Polycarbonate Body, fast disinfection process, protects UV and IR rays and improve light penetration. Aerodynamic design for ventilation of heat. Body is made of high quality fire resistant polymer



Memory Function



Endo Mode



Battery Backup: Optional

## Star 63

Intensity	1,60,000 Lux
Shadowless Feature	Yes
Numbers of LED's	63
Intensity Control	10~100 Digital Capacitive Touch (Fully Remote Controlled)
Fitting Options Available	Ceiling 180 Degree, Ceiling 360 Degree, Mobile Stand, Wall Mounted
Spot Diameter	100-150 mm
Depth of Light Penetration	150 mm
Colour Temperature	3800 – 4800 K
Colour Rendering Index	95 Ra
Focusing	Adjustable
LED Average life	>50,000 Hrs.
Diameter of Light	700 mm
Power Supply	220V/50Hz AC
Power Consumption	100 Watt

## Features

Shadowless Light, Memory Function, High quality sterilizable handle, With protective fire resistant Polycarbonate Body, fast disinfection process, protects UV and IR rays and improve light penetration. Aerodynamic design for ventilation of heat. Body is made of high quality fire resistant polymer



Memory Function



Endo Mode



Battery Backup: Optional