

# LED EXAMINATION LIGHTS

Examination Lights are used for Minor Operations, Stitches and Examination Purpose.

Our diverse range of Examination Lights are very user-friendly as they are easy-to-use, easy to install and help in achieving the perfect viewing angle as they are easy to position using Flexible Gooseneck Pipe.



**Examination lights for minor surgery purpose  
Detachable Mobile Stand**

## MAGNA

Intensity	60,000 Lux
Special Feature	Inbuilt magnifying glass
Numbers of LEDs	8 (Warm + Cool White)
Intensity Control	10-100 Digital Capacitive Touch
Fitting Options Available	Mobile Mounted
Spot Diameter	100-150 mm
Depth of Light Penetration	100 mm
Colour Temperature	4800 K
Colour Rendering Index	93 Ra
Focusing Preset	at 1 Metre
LED Average life	>50,000 Hrs.
Diameter of Light	160 mm
Power Supply	220V/50Hz AC
Power Consumption	15 Watt



## Features

Magnification glass inbuilt for micro-zooming during surgeries, stitches. Perfect solution for minor surgeries, hair transplant, dermatologists, gynaecologists etc.



Battery Backup: Optional

## Venus 1

Intensity	50,000 Lux
Numbers of LEDs	24 (Warm + Cool White)
Intensity Control	10-100 Digital Capacitive Touch
Fitting Options Available	Mobile Stand
Spot Diameter	100-150 mm
Depth of Light Penetration	100 mm
Colour Temperature	4800K
Colour Rendering Index	93 Ra
Focusing Preset	at 1 Metre
LED Average life	>50,000 Hrs.
Diameter of Light	150 mm
Power Supply	220V/50Hz AC
Power Consumption	20 Watt



## Features

Unbreakable Polycarbonate Body | Touch Controller



Battery Backup: Optional

## Star 21

Intensity	50,000 Lux
Spot Diameter	100-150 mm
Depth Diameter	10 -25 cm
Color Temp.	4000k
Colour Rendering Index	93 Ra
Focusing	Manual
LED Average life	>50000 Hrs.
Numbers of LED's	21
Diameter of Light	200 mm
Intensity Control	10-100 Digital Touch
Power Supply	220V/50Hz AC
Power Consumption	15 Watt
Special Function	Feather Touch Digital Controller



## Features

Unbreakable Polycarbonate Body | Touch Controller



Battery Backup: Optional

## PLUTO 07

Intensity	25,000 Lux
Numbers of LEDs	7
Intensity Control	On/off Button
Fitting Options Available	Mobile mounted
Colour Temperature	4500 K
Colour Rendering Index	93 Ra
Focusing Preset	0.5 Metre
LED Average life	>50,000 Hrs.
Diameter of Light	150mm
Power Supply	220V/50Hz AC
Power Consumption	10 Watt



## Features

ON/OFF Button | Unbreakable Polycarbonate Body



Battery Backup: Optional