

अनुबंध | Contract



अनुबंध क्रमांक | Contract No: GEMC-511687715482158

अनुबंध तिथि | Generated Date : 10-Jun-2025

संगठन विवरण Organisation Details		खरीदार विवरण Buyer Details				
प्ररूप Type :	Central Government	पद Designation :	Assistant Professor Pediatrics Department			
मंत्रालय Ministry :	Ministry of Health and Family Welfare	संपर्क नंबर Contact No. :	-			
विभाग Department :	Department of Health and Family Welfare	ईमेल आईडी Email ID :	buycon88.aiimsa.hp@gembuyer.in			
संगठन का नाम Organisation Name :	All India Institute of Medical Sciences (AIIMS)	जीएसटीआईएन GSTIN :	-			
कार्यालय क्षेत्र Office Zone:	bilaspur himachal pradesh	पता Address :	aiims ,kothipura,changar palasia, Bilaspur, HIMACHAL PRADESH-174001, India			
वित्तीय स्वीकृति विवरण Financial Approval Detail		भुगतान प्राधिकरण विवरण Paying Authority Details				
आईएफडी सहमति IFD Concurrence :	Yes	Role:	PAO			
प्रशासनिक अनुमोदन का पदनाम Designation of Administrative Approval:	Executive Director , AIIMS, Bilaspur , H.P	भुगतान का तरीका Payment Mode:	Offline			
वित्तीय अनुमोदन का पदनाम Designation of Financial Approval :	Executive Director , AIIMS, Bilaspur , H.P	पद Designation :	Accounts Officer			
		ईमेल आईडी Email ID :	pay5.aiimsa.hp@gembuyer.in			
		जीएसटीआईएन GSTIN :	02PTLA18254F1DX			
		पता Address:	aiims ,kothipura,changar palasia, Bilaspur, HIMACHAL PRADESH-174001, India			
विक्रेता विवरण Seller Details						
जेम विक्रेता आईडी GeM Seller ID :	24D1180000472836					
कंपनी का नाम Company Name :	Vandana Surgicals					
संपर्क नंबर Contact No. :	09999049518					
ईमेल आईडी Email ID :	Keshavsinghal25@gmail.com					
पता Address :	BH Block,Vandana Surgicals,BH 14, East, Shalimar Bagh, New Delhi 110088,Jhulelal Mandir Marg, North West delhi, DELHI-110088, -					
एमएसएमई पंजीकरण संख्या MSME Registration number :	UDYAM-DL-04-0002130					
जीएसटीआईएन GSTIN:	07AAWPS1302H1ZK (R) , (M)					
एमएसई सामाजिक श्रेणी MSE Social Category :	General					
एमएसई लिंग श्रेणी MSE Gender :	Female					
*जिसके नाम के पक्ष में GST/TAX इनवॉइस पेश किया जाएगा GST / Tax invoice to be raised in the name of - Consignee						
वितरण निर्देश Delivery Instructions : NA						
उत्पाद विवरण Product Details						
#	आइटम विवरण Item Description	आइटम विवरण Ordered Quantity	इकाई Unit	इकाई मूल्य (INR) Unit Price (INR)	कर विभाजन (INR) Tax Bifurcation (INR)	मूल्य (INR में सभी शुल्क और कर सहित) Price (Inclusive of all Duties and Taxes in INR)
1	उत्पाद का नाम Product Name : Ventek India Examination Light for Medical Use ब्रांड Brand : Ventek India ब्रांड प्रकार Brand Type : Registered Brand कैटलॉग की स्थिति Catalogue Status: OEM verified catalogue कैसे बेचा जा रहा है Selling As : OEM श्रेणी का नाम और चतुर्थांश Category Name & Quadrant : Examination Light for Medical Use (V2) (Q2) मॉडल Model: Star 21 एचएसएन कोड HSN Code: HSN not specified by seller	1	pieces		NA	
कुल ऑर्डर मूल्य Total Order Value (in INR)						
परोक्षिती विवरण Consignee Detail						
क्र.सं. S.No	परोक्षिती Consignee	वस्तु Item	लॉट नंबर Lot No.	मात्रा Quantity	दिनांक के बाद डिलीवरी शुरू करना है Delivery Start After	वितरण पूरा कब तक करना है Delivery To Be Completed By
1	पद Designation : Assistant Professor Pediatrics Department ईमेल आईडी Email ID : buycon88.aiimsa.hp@gembuyer.in संपर्क Contact : - जीएसटीआईएन GSTIN : -	Ventek India Examination Light for	-	1	10-Jun-2025	25-Jun-2025

पता Address : aiims ,kothipura,changar palasia, Bilaspur, HIMACHAL PRADESH-174001, India		Medical Use				
Product Specification for Ventek India Examination Light for Medical Use						
विनिर्देश Specification		उप-विनिर्देश Sub-Spec			मूल्य Value	
GENERAL	Product Name			Examination Light for Medical Use		
	Purpose			Examination light provides bright, cool, and adjustable illumination to enhance visibility during medical examinations and minor procedures		
LIGHT SOURCE	Light Type			LED		
	Correlated Colour Temperature (in Kelvin)			3000 to 4500 K		
	Colour Rendering Index (CRI)			80 to 100		
	Working Distance (in mm)			100 to 400 mm		
	Minimum Light Intensity at working distance (in Lux)			50000		
	Spot size at working distance (in mm)			60 to 150 mm		
	Minimum LED Life (in Hours)			50000		
POWER & CONTROLS	Power Supply			220 V +/- 10%, 50 Hz		
	Availability of continuous dimming levels			Yes		
	Touchless On/Off control			Yes		
	Length of the Power Chord (in meters)			>=1		
	Power Consumption (in watts)			7 to 150 W		
MOUNTING & ACCESORIES	Mounting options			Mobile Stand		
	Wheel locking provision			Yes		
	Provision for extension (Goose Neck)			Yes		
CERTIFICATIONS	Compliance to Medical Device Rules (MDR) 2017 as amended till date			Yes		
	Availability of valid Medical Device license for the product issued from the competent authority defined under Drugs and Cosmetic Act 1940 and Rules made there under as amended till date			Yes		
	Valid Medical Device License Number			MFG/MD/2023/000867		
	Certification for manufacturing unit			ISO:13485 (Latest)		
	Availability of Test Report for each supplied batch/product as per Medical Device Rule (MDR) 2017 as amended till date			Yes		
	Submission of all necessary certifications, licenses and test reports to the buyer at the time of bid submission or along with supplies as per buyer requirement			Yes		
	Electrical Safety Compliance Standard			IEC 60601 or BIS Equivalent		
WARRANTY	Warranty in years (Option of comprehensive warranty is available through bidding only, which if opted will supersede normal warranty in the catalogue)			2 year		
ADDITIONAL REQUIREMENTS	Additional Requirements			NA		
ईपीबीजी विवरण ePBG Detail						
NA						
नियम और शर्तें Terms and Conditions						
1. Special terms and conditions- Version:1 effective from 30-08-2024						
1.1	1. All Provisions of Drugs and Cosmetics Act, 1940 and Rules (including Medical Device Rule 2017) made there under as amended till date will always be applicable. This will include all notifications issued by Central Drugs Standard Control Organization (CDSCO), Ministry of Health & Family Welfare (MoHFW) and Department of Pharmaceuticals (DOP), Ministry of Chemicals & Fertilizers time to time in this regard. 2. The sellers are registered on GeM based on self declaration of valid Medical Device License, product certification, test reports etc. However, buyers must check and validate the details at their end for all applicable licenses and certifications e.g., validity and authenticity/genuineness of Medical Device license, product certification, manufacturer certification/licenses, test reports etc. 3. In case of authorized resellers/distributors, it will be the legal & regulatory liability of the manufacturer to ensure that their resellers/distributors are operating in compliance with all relevant laws and regulations and are properly licensed to sell the manufacturer's products, including verifying the validity and authenticity of Medical Device license held by them. 4. The price offered by the seller/bidder shall not, in any case exceed the DPCO/NPPA controlled price or price fixed by State Government, if any. The seller must reduce					

the prices if there is any reduction in DPCO/NPPA ceiling price or price fixed by State Government, if any.

5. Any other Terms and Conditions which is not included or at variance with the conditions specified in STC/GTC, may be added by the buyer through Additional Terms and Conditions (ATC) in the bid to ensure items are procured from authentic/validated source with appropriate and applicable quality. The above terms and conditions are in reverse order of precedence i.e. ATC shall supersede specific STC which shall supersede General Terms and Conditions (GTC), whenever there are any conflicting provisions.
6. Comprehensive warranty: Comprehensive warranty shall include preventive maintenance including calibration as per technical/ service /operational manual of the manufacturer, service charges and spares. During the warranty period commencing from date of the successful completion of warranty period, Service personnel shall visit each consignee site as recommended in the manufacturer's technical/ service /operational manual, at least once in six months. warranty shall not be including the consumables. Further there will be 98% uptime warranty during warranty period on 24 (hrs) X 7 (days) X 365 (days) basis, with penalty, to extend warranty period by double the downtime period.
7. Service centres: Details of Service outlets in India to render services for equipment to be furnished to buyer/consignees with complete address, telephone numbers, e mails etc at time of making the supplies. It shall be the responsibility of seller to ensure that authorized service centres are available to cater to the areas where supplies are made within reasonable distance from where the service calls can be handled. Details of toll-free numbers for service call and online registration of service requests also to be provided buyer/consignee at the time of supplies.
8. Source of supply: It shall be responsibility of seller to provide Documents regarding source of equipments such as copy of Performa invoice or any other documents to establish that the products supplied are manufactured by OEM indicated and sourced from them.
9. Packing and Marking: Medical equipments being very delicate and sensitive packing for the goods should be strong and durable enough to withstand transit including transshipment (if any), rough handling, open storage etc. without any damage, deterioration etc. .The size, weights and volumes of the packing cases, remoteness of the final destination of the goods, availability or otherwise of transport and handling facilities at all points during transit up to final destination,. Quality of packing, the manner of marking within & outside the packages and provision of accompanying documentation shall take into consideration the type of medical equipments being supplied. The accessories shall be suitably labelled and packed. Each of the package shall be marked on three sides with indelible paint of proper quality: indicating contract number and date, brief description of goods including quantity, Packing list reference number, country of origin of goods and any other relevant details.
10. Spare Parts: Seller shall provide materials, information etc. pertaining to spare parts manufactured and supplied by the OEM. It shall be ensured that the required spares are available for purchase at least for 10 years from date of supplies. In case due to any reasons the production of the spare parts is discontinued sufficient advance notice should be given to the buyer/consignee before such discontinuation to provide adequate time to purchase the required spare parts etc. Further, OEM and their service centres/dealers shall carry sufficient inventories to assure ex-stock supply of consumables and spares for the equipments so that the same are available. OEM or reseller shall always accord most favoured client status to the buyer/consignee and shall give the most competitive price for spares and consumables of its machines/equipments supplied.
11. Installation, Training, Manuals: Seller shall be responsible to carry out Installation & commissioning, Supervision and Demonstration of the goods. They shall provide required jigs and tools for assembly, minor civil works for the completion of the installation and Training of Consignee's representatives for operating and maintaining the equipment and supplying required number of operation & maintenance manual for the goods. In case the category parameters are specifying any requirements regarding the installations, training and manuals the same shall also be applicable.
12. Electrical safety checking: Sellers are required to make sure that they furnish the list of equipments for carrying out routine and preventive maintenance to buyer/consignee .They should make sure to periodically check the electrical safety aspects as per BIS Safety Standards or equivalent .In case they do not have required equipment for such testing should ensure that the equipments checked for electrical safety compliance through labs with facilities for such checking during every preventive maintenance call.
13. Software: All software updates should be provided free of cost during warranty period.

2. General Terms and Conditions-

- 2.1 This contract is governed by the [General Terms and Conditions](#), conditions stipulated to this Product/Service as provided in the Marketplace.
- 2.2 This Contract between the Seller and the Buyer, is for the supply of the Goods and/ or Services, detailed in the schedule above, in accordance with the General Terms and Conditions (GTC) unless otherwise superseded by Goods / Services specific Special Terms and Conditions (STC) and/ or BID/Reverse Auction Additional Terms and Conditions (ATC), as applicable
- 2.3 All GeM Sellers / Service Providers are mandated to ensure compliance with all the applicable laws / acts / rules including but not limited to all Labour Laws such as The Minimum Wages Act, 1948, The Payment of Wages Act, 1936, The Payment of Bonus Act, 1965, The Equal Remuneration Act, 1976, The Payment of Gratuity Act, 1972 etc. Any non-compliance will be treated as breach of contract and Buyer may take suitable actions as per GeM Contract.

नोट: यह सिस्टम जनरेटेड फाइल है। कोई हस्ताक्षर की आवश्यकता नहीं है। इस दस्तावेज़ का प्रिंट आउट भुगतान/लेनदेन उद्देश्य के लिए मान्य नहीं है।

Note: This is system generated file. No signature is required. Print out of this document is not valid for payment/ transaction purpose.